MACHEREY-NAGEL

Rapid Tests for Schools



- Economical
- High educational value
- Easy to use





MN Dip & Read Tests for Schools

Test strips motivate kids to discover their environment in many different ways. They can easily check the water quality in the local river, the ozone concentration in the summer air, the vitamin C content of fruits etc. The easy Dip & Read procedure of test strips safely provides reliable results in a few seconds, making MN products the perfect choice for an exiting but safe educational environment.



Product information

Product	Application	Gradation	REF
Ozone Test strip	 Ozone is formed on sunny days from oxygen and nitrogen oxides 	< 90 · 90 - 150 · 150 - 210 · > 210 μg/m³ O ₃	90736
	 Easy test to measure the local ozone concentration 		
	Provides		
AQUADUR® 3-25	■ Easy determination of water hardness	< 3 · > 5 · > 10 · > 15 · > 20 · > 25 °d	91201
	■ 5 different pads easy read off		
	 Provides valuable information about the local water 		
QUANTOFIX® Nitrate / Nitrite	Fast and easy detection of nitrate and nitrite	$0 \cdot 10 \cdot 25 \cdot 50 \cdot 100 \cdot 250 \cdot 500 \text{ mg/l NO}_3^-$ $0 \cdot 1 \cdot 5 \cdot 10 \cdot 20 \cdot 40 \cdot 80 \text{ mg/l NO}_2^-$	91313
	 Check water quality in the local river 		
	 Detects nitrate in vegetables and potatoes 		
pH-Fix 0 – 14	 Easy detection of pH value in surface water, rivers, fruit juice etc. 	pH 0 · 1 · 2 · 3 · 4 · 5 · 6 · 7 · 8 · 9 · 10 · 11 · 12 · 13 · 14	92110
	 Brilliant colour charts for clear results 		
	 Plastic strip prevents skin contact with the sample 		
Others			
Litmus paper, neutral	Basic determination of acids and bases		91107
Potassium iodide starch paper	Detection of oxidisers like free chlorine and nitrite		90754
QUANTOFIX® Ascorbic acid	For vitamin C in juice or fruits		91314
Oil test paper	Oil detection in water		90760
Sulfide test paper	Detection of hydrogen sulfide		90761
Medi-Test Glucose	Specific test for glucose		93024
Watesmo	Easy detection of water		90609

Your local distributor





Follow us



www.mn-net.com

HEREY-NAGEL



MACHEREY-NAGEL GmbH & Co. KG Valencienner Str. 11 52355 Düren · Germany

DE Tel.: +49 24 21 969-0

info@mn-net.com CH Tel.: +41 62 388 55 00 sales-ch@mn-net.com FR Tel.: +33 388 68 22 68 sales-fr@mn-net.com US Tel.: +1 888 321 62 24 sales-us@mn-net.com

