

## EZ-CFU™ ONE STEP

### FOR GROWTH PROMOTION TESTING

- Quantitative
- Refrigerated storage
- Eight hour stability
- No dilutions required
- Peel-off Certificate of Assay
- Online Certificate of Analysis
- Traceable to reference culture
- Product warranty



## EZ-CFU™ ONE STEP PRODUCT DETAILS

The Growth Promotion Test (GPT) is an important function in the United States (USP), European (Ph. Eur.) and Japanese (JP) Pharmacopeias. EZ-CFU™ One Step quantitative microorganism preparations are used to perform Growth Promotion Tests of culture media with ease. No dilutions are required. EZ-CFU™ One Step is designed to deliver less than 100 CFU with each 0.1 mL inoculum. Each 2 mL vial of hydrated suspension offers 19 inocula. EZ-CFU™ One Step can be used for up to eight hours after hydration if the suspension is refrigerated between uses.\* For added convenience, a peel-off Certificate of Assay is included for quality control documentation.

\*Catalog numbers 0318Z and 0320Z are an exception to the eight hour stability and must be used within 30 minutes of hydration.

Distribué par :

**LABORATOIRES HUMEAU**

Z. A. de Gesvrine - 4 rue Képler - B. P. 4125 - 44241 La Chapelle-sur-Erdre Cedex - France  
t. : +33 (0)2 40 93 53 53 - f. : +33 (0)2 40 93 41 00 - e. : info@humeau.com