

EZ-CFU™

FOR GROWTH PROMOTION TESTING

- Quantitative
- Refrigerated storage
- Peel-off Certificate of Assay
- Online Certificate of Analysis
- Traceable to reference culture
- For high volume users
- Product warranty



EZ-CFU™ PRODUCT DETAILS

EZ-CFU™ is used to perform Growth Promotion Tests of culture media with ease. EZ-CFU™ lyophilized microorganism preparations deliver less than 100 CFU per 0.1 mL when processed as directed. A one log dilution in pH 7.2 phosphate buffer is required. Following the dilution step, EZ-CFU™ provides over 90 inocula using the same suspension. For added convenience, a peel-off Certificate of Assay is included for quality control documentation.

Distribué par :

LABORATOIRES HUMEAU

Z. A. de Gesvrine - 4 rue Képler - B. P. 4125 - 44241 La Chapelle-sur-Erdre Cedex - France
t. : +33 (0)2 40 93 53 53 - f. : +33 (0)2 40 93 41 00 - e. : info@humeau.com



w w w . h u m e a u . c o m