



KWIK-STIK™ PLUS

TWO PASSAGE MICROORGANISMS

- Allows more subculturing
- Refrigerated storage
- Peel-off Certificate of Assay
- Online Certificate of Analysis
- Traceable to reference culture
- Product warranty

KWIK-STIK™ PLUS PRODUCT DETAILS

KWIK-STIK™ Plus microorganisms are designed to meet the needs of laboratories requiring QC microorganisms that are a limited number of passages from the reference culture. Each KWIK-STIK™ Plus microorganism preparation is just two passages from the reference culture, allowing for more subcultures. Conveniently packaged in a KWIK-STIK™ device featuring a lyophilized microorganism pellet, inoculating swab and reservoir of hydrating fluid, the KWIK-STIK™ Plus is ready to use and affordable. For added convenience, a peel-off Certificate of Assay is included for quality control documentation. KWIK-STIK™ Plus includes five STIKS of a single strain.

KWIKSTIK™
plus

Distribué par :

LABORATOIRES HUMEAU

Z. A. de Gesvrine - 4 rue Képler - B. P. 4125 - 44241 La Chapelle-sur-Erdre Cedex - France
t. : +33 (0)2 40 93 53 53 - f. : +33 (0)2 40 93 41 00 - e. : info@humeau.com

w w w . h u m e a u . c o m