

KWIK-STIK™

FOR CULTURE PURPOSES AND QC APPLICATIONS

- Economical
- Refrigerated storage
- Online Certificate of Analysis
- Traceable to reference culture
- Product warranty



KWIK-STIK™ PRODUCT DETAILS

The KWIK-STIK™ is well-known for its clever design and accurate, repeatable results. Each KWIK-STIK™ device features a single microorganism strain in a lyophilized pellet, a reservoir of hydrating fluid, and inoculating swab. For added convenience, a peel-off identification label is included for easy documentation. The KWIK-STIK™ is available in packs of 2 or 10.

KWIKSTIK™

Distribué par :

LABORATOIRES HUMEAU

Z. A. de Gesvrine - 4 rue Képler - B. P. 4125 - 44241 La Chapelle-sur-Erdre Cedex - France
t. : +33 (0)2 40 93 53 53 - f. : +33 (0)2 40 93 41 00 - e. : info@humeau.com



w w w . h u m e a u . c o m