

KWIK-STIK™

FOR CULTURE PURPOSES AND QC APPLICATIONS

- Economical
- Refrigerated storage
- Online Certificate of Analysis
- Traceable to reference culture
- Product warranty



KWIK-STIK™ PRODUCT DETAILS

The KWIK-STIK™ is well-known for its clever design and accurate, repeatable results. Each KWIK-STIK™ device features a single microorganism strain in a lyophilized pellet, a reservoir of hydrating fluid, and inoculating swab. For added convenience, a peel-off identification label is included for easy documentation. The KWIK-STIK™ is available in packs of 2 or 10.

