
1x20 SACHET 25ml EC 1413 $\mu\text{S}/\text{cm}$ $\pm 1\%/25^\circ\text{C}$, with Certificate of Analysis

Réf.: 203386



DATA SHEET

Features

Accuracy	$\pm 1\%$
Certificate	Product Reports
Colour	Colourless
Composition	Potassium chloride in aqueous solution
Pack	20 sachets of 25 ml
Description	Box of 20 sachets of 25 ml, with certificate of analysis
Image 2	https://app.onpage.it/api/storage/500c8ce02e5c84b5b7bd34fdb79bfaf76ee53
Name	EC 1413 $\mu\text{S}/\text{cm}$
Nominal value	1413 $\mu\text{S}/\text{cm}$ (25 $^\circ\text{C}$)
Volume	25 ml

DS_51102053_rev._18-10-2023



Distribué par :

Z.A de Gesvrine - 4 rue Képler - B.P.4125
44241 La Chapelle-sur-Erdre Cedex - France
t. : +33 (0)2 40 93 53 53 | f. : +33 (0)2 40 93 41 00
commercial@humeau.com