
1X500ml XS Solution pH 10.01 ±0.02/25°C (colourless), with certificate traceable to NIST

Réf.: 203247



DATA SHEET

Features

Accuracy	± 0.02
Certificate	Yes
Colour	Colourless
Composition	Sodium carbonate, sodium bicarbonate
Pack	500 ml with integrated calibration compartment
Description	±0.02 / 25°C - Colourless
Name	pH 10.01
Nominal value	pH 10.01 @ 25 °C
Volume	500 ml

DS_51100163_rev._18-10-2023



Distribué par :

Z.A de Gesvrine - 4 rue Képler - B.P.4125
44241 La Chapelle-sur-Erdre Cedex - France
t. : +33 (0)2 40 93 53 53 | f. : +33 (0)2 40 93 41 00
commercial@humeau.com