

FP4

The **FP4** is the **portable flocculator** that can be connected to a battery or a car cigarette lighter.

The anti-skid base ensures **stability** whilst the highly versatile stainless steel stirring shafts offer height adjustment and are fitted with a self-locking device with clutch.

The FP4 is particularly suitable for carrying out on-site Jar Tests.

The **carrying case**, an optional accessory, is ideal for storing or transporting the unit.

UK, AU and USA adapter plugs are available on request.



INSTRUMENT	POWER SUPPLY	CODE No
FP4	100-240 V / 50-60Hz	F105A0117

GENERAL FEATURES

Construction material:	Epoxy painted metal structure
Disconnectable central light	
Stainless steel stirring rods:	adjustable in height by a self blocking chuck
DC gear motor	
Monoselector:	same speed for each rod
Number of stirring rods:	4
Power:	6 W
Power supply:	100 - 240 V / 50-60 Hz
Weight:	4,8 Kg (10.6 lb)
Dimensions (WxHxD) :	250x320x250 mm (9.8x12.6x9.8 in)

PERFORMANCE

Timer:	from 0 to 30 minutes (or continuous)
Electronic speed control:	20 - 40 - 50 - 100 - 200 rpm

Distribué par :

Z.A de Gesvrine - 4 rue Képler - B.P.4125
44241 La Chapelle-sur-Erdre Cedex - France
t. : +33 (0)2 40 93 53 53 | f. : +33 (0)2 40 93 41 00
commercial@humeau.com

