

GANTS NITRILE BLEU NON POUDRÉS



General information:

Type:	disposable gloves
Description:	All information is printed on the dispenser box.
Shape:	ambidextrous, straight fingers
Material:	Nitrile Butadiene Rubber (NBR)
Color:	blue (recommended for HACCP Standard)
Inner:	online chlorinated, powder free
Outer:	no treatment
Cuff/Surface:	rolled cuff, fully textured surface
Length:	Ø 240 mm
Thickness (double):	Finger min. 0.25 mm, Palm min. 0.25 mm, Cuff min. 0.16 mm
Force at break (EN 455-2):	Ø 6,0 N (during shelf life)
Tensile strength/elongation:	14 MPa / 500 % (after aging: 14 MPa / 400 %)
Residual powder:	≤ 2 mg / gloves

Available sizes: S (6 - 7), M (7 - 8), L (8 - 9), XL (9 - 10)

Storage:

Storage life:	3 years (expiry date is printed on the dispenser box).
Storage conditions:	cool & dry (recommended: 10 - 30 °C at 60 - 80% humidity), dust free, avoid UV radiation, avoid contact with oxidizing agents and ozone

Directives and Certification:

Directives:	89/686/EEC Personal Protective Equipment, CE category III
Standards:	EN 420, EN 455, EN 374-2 (specific AQL value 1.5), EN 374-3 (Class 6, Sodium hydroxide 40%).
Food regulation:	1935/2004

Warning:

Gloves contain the following allergenic substances: Dithiocarbamate

This product contains allergenic substances which can cause allergic reactions, including anaphylactic responses.