



Product Specifications

CL00.2301

Hydrogen peroxide 35 weight % solution (stab.) p.

Waterstofperoxide 35 gew. % oplossing (stab.) z.
Hydrogène peroxyde 35 poids % solution (stab.) p.
Wasserstoffperoxid 35 w/w % Lösung (stab.) Rn.
Hidrógeno peróxido 35 p./v % solución (stab.) p.
Perossido di idrogeno 35 p/p % soluzione (stab.) p.

For laboratory use
35+% H₂O₂ (phosphate stabilised)

Mol weight: 34.01 g/mol
Density: 1,13 g/ml
HS Nr: 28470000
CAS Nr: 7722-84-1
UN: 2014
ADR: 5.1 (8),II



Specifications:

Assay	>35%
Non Volatiles	<0.05%
Residue after Ignition	<0.03%
Free Acid	<0.02%
Heavy Metals as Lead (Pb)	<0.0002%
Phosphate (PO ₄)	<0.05%