

---

# M2-A ARGOLab MAGNETICAL STIRRER- Analogic Control

Item Code: 22002013



## DATA SHEET

---

### Features

Power supply	220 V
Dimension	280 x 160 x 65 mm
Frequency	50 / 60 Hz
Warranty	2 anni
IP rating	42 IP
Temperature indication	Analogue
Speed indication	Analogue
Max lenght magnetic bar	80 mm
Motor	Brushless DC (Direct Current) motor
Name	M2-A
Weight	2.8 kg
Working plate - dimensions	Ø 135 mm
Working plate - material	Stainless Steel
Heating power	500 W
Over Temperature Protection	420
Working temperature	5 ... 40 °C
Heating Temperature Range	Ta ... 340 °C
Working humidity	80 %RH
Stirring speed	0 ... 1500 rpm
Max stirring quantity	20 L

---

DS\_22002013\_rev. \_18-01-2024