

# M2-D PRO ARGOLab

## MAGNETICAL

## STIRRER- Digital Control

Item Code: *22002063*

## DATA SHEET

### Features

Accuracy with external probe	± 0,2 °C
"Hot" warning	50 °C
Connection	RS 232
Heating accuracy control	±1 °C (<100 °C) - ±1 % (>100 °C)
Dimension	280 x 160 x 65 mm mm
Max anchor size	80 mm
Frequency	50 / 60 Hz
Warranty	2 anni
IP rating	42 IP
Image 2	<a href="https://app.onpage.it/api/storage/a6e06e832c6b60c42c0e6e7a6f61d2c476ee0">https://app.onpage.it/api/storage/a6e06e832c6b60c42c0e6e7a6f61d2c476ee0</a>
Temperature indication	Digital
Speed indication	LCD Display
Max lenght magnetic bar	80 mm
Motor	Brushless DC (Direct Current) motor
Name	M2-D Pro
Weight	2.8 kg
Working plate - heating area	135 mm
Working plate - dimensions	Ø 135 mm
Working plate - material	Ceramic-coated stainless steel
Absorbed power	550 W
Heating power	500 W
Over Temperature Protection	50
External temperature probe	PT 1000
Working temperature	5 ... 40 °C
Heating Temperature Range	Ta ... 340 °C °C
Working humidity	80 %RH
Stirring speed	0 ... 1500 rpm
Max stirring quantity	20 L