

# M3-D Agit. ARG0lab Magnet.Stirrer with heating ceramic plate

Item Code: 22003013



## DATA SHEET

### Features

Accuracy with external probe	± 0,2 °C
Power supply	220 V
Connection	RS 232
Dimension	215 x 360 x 112 mm
Max anchor size	80 mm
Frequency	50 / 60 Hz
Warranty	2 anni
IP rating	21 IP
Image 2	<a href="https://app.onpage.it/api/storage/317f81c74c06b76116d6df5605fb298a320def!">https://app.onpage.it/api/storage/317f81c74c06b76116d6df5605fb298a320def!</a>
Temperature indication	LED Display
Speed indication	Digital
Max lenght magnetic bar	80 mm
Motor	Shielded rod motor
Name	M3-D
Weight	4.5 kg
Working plate - dimensions	Ø 135 mm
Working plate - material	Glass ceramic
Absorbed power	1030 W
Heating power	1000 W
Over Temperature Protection	580
External temperature probe	PT 1000
Working temperature	Ta ... 550 °C
Heating Temperature Range	Ta ... 550 °C °C
Working humidity	80 %RH
Stirring speed	0 ... 1500 rpm
Max stirring quantity	10 L