
MS-D ARG0lab digital magnetical stirrer

Item Code: 22002003



DATA SHEET

Features

Accuracy (stirrer speed)	±1 rpm
Power supply	220 V
Dimension	260 x 150 x 80 mm mm
Frequency	50 / 60 Hz
Warranty	2 anni
IP rating	42 IP
Speed indication	Digital via LED display
Max lenght magnetic bar	25 mm
Name	MS-D
Weight	1.8 kg
Working plate - dimensions	135 mm
Working plate - material	Nylon+GF
Absorbed power	15 W
Motor input power	5 KW
Motor output power	3 KW
Working humidity	80 %RH
Stirring speed	200 - 1500 rpm
Max stirring quantity	3 L

DS_22002003_rev. _18-01-2024



Distribué par :

Z.A de Gesvrine - 4 rue Képler - B.P.4125
44241 La Chapelle-sur-Erdre Cedex - France
t. : +33 (0)2 40 93 53 53 | f. : +33 (0)2 40 93 41 00
commercial@humeau.com