

Product Specifications

CL00.0165

Amylalcohol-(iso) v.p.

Alcole amílico-(iso) e.p.
Amylalcohol-(iso) Rst.
Alcool amylique-(iso) t.p.
Alcohol Amílico-(iso) t.p.
Amylalcohol-(iso) z.z.

Mol weight: 88.15 g/mol
Density: 0,81 g/ml
HS Nr: 29051900
CAS Nr: 123-51-3
UN: 1105
ADR: 3,III



For laboratory use, Fat det. acc. Gerber, NF V 04-210
98+% C₅H₁₂O (±10% 2-Methyl-1-butanol and ± 90% 3-methyl-1-butanol)

Specifications:

Assay	>98%
Water	<0.3%
Density	= 0.808 - 0818 g/ml
Distillation Range	= 128-132°C
2-Methyl-1-butanol	= 8-10%
3-Methyl-1-butanol	= 90-92%
Colour	< 10 APHA