

CL00.2722

CHROMALYT

Ammonium acetate, eluent additive for LC-MS

Ammonio acetato, LC-MS
Ammoniumacetat, LC-MS
Amonio acetato, LC-MS
Ammonium acétate, LC-MS
Ammoniumacetaat, LC-MS

Mol weight: 77.08 g/mol
Density: 1,17 g/cm³
HS Nr: 29152900
CAS Nr: 631-61-8
UN:
ADR:

For laboratory use, eluent additive for LC-MS
99+% C₂H₇NO₂

Specifications:

Assay	>99%
Ammonium Acetate	>99%
Aluminium (Al)	<0.0001%
Barium (Ba)	<0.0001%
Cadmium (Cd)	<0.0001%
Calcium (Ca)	<0.001%
Chromium (Cr)	<0.0001%
Cobalt (Co)	<0.0001%
Copper (Cu)	<0.0001%
Iron (Fe)	<0.0001%
Lead (Pb)	<0.0001%
Lithium (Li)	<0.0001%
Magnesium (Mg)	<0.0001%
Manganese (Mn)	<0.0001%
Molybdenum (Mo)	<0.0001%
Nickel (Ni)	<0.0001%
Potassium (K)	<0.005%
Sodium (Na)	<0.005%
Strontium (Sr)	<0.0001%
Zinc (Zn)	<0.0001%
Satisfies LC-MS applications	