

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Anti-Bovine Albumin antibody produced in rabbit - whole antiserum

Product Number:

B1520

TEST	Specification
------	---------------

Appearance (Form)

Liquid

Appearance (Turbidity)

Clear to Slightly Hazy

Appearance (Color)

Very Light Brown to Very Dark Brown and
Dark Brown-Red to Very Dark Brown-Red

Protein by Biuret

40 - 100 mg/ml

Antibody Concentration(Method)

Quantitative Precipitin Assay

Antibody Concentration

3.0 - 5.0 mg/ml

ELISA - Titer

1:200,000-1:300,000

Recommended Retest Period

5 Years

Specification: PRE.1.ZQ5.10000016050

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.