

Technical Data Sheet

Doc No – Ab/TDS/Micro-83-I

Rev No: 07 Date 05/2024-25

TECHNICAL DATA SHEET

Abdos Labtech Pvt. Ltd.

Petri Dish, PS - Aseptic

Valid for Cat. No. P10901, P10902, P10903, P10904, P10931, P10932, P10933, P10934

SRN (Single Registration Number): IN-MF-000022272

1. Description/Specification

1.1 Description: Petri Dish- Aseptic, DNase, RNase, Human DNA, ATP , Pyrogen Free, CE & IVD certified, ISO24998

1.2 Design: Shallow dish design

1.3 Performance Specifications:

Code	Diameter (mm)	Height (mm)	Compartment	Vented	Packing	No. Per Case
P10901	90	15.80	No	3 Vent	Sleeve of 20	480
P10902	90	15.80	No	3 Vent	Individual Wrapped	380
P10903	90	15.80	No	0 Vent	Sleeve of 20	480
P10904	90	15.80	No	0 Vent	Individual Wrapped	380
P10931	60	14.20	No	0 Vent	Sleeve of 20	480
P10932	60	14.20	No	3 Vent	Sleeve of 20	480
P10933	140	17.50	No	0 Vent	Sleeve of 5	180
P10934	140	17.50	No	3 Vent	Sleeve of 5	180

1.4 Material:

1.4.1 Petri Dish made up of Hydrophobic Polystyrene (PS), Confirm with US FDA 21CFR

1.4.2 Complies with the requirement of European Pharmacopoeia, Plastic Containers and closures for pharmaceutical use.

1.4.3 The product has passed the USP testing including class VI tests, Materials has successfully passed the biological tests according to ISO 10993 external communication devices for indirect blood contact for Prolonged period.

1.4.4 Free from Heavy metals, RoHS Free according to meets relevant requirements of EU regulations 1935/2004/CE, 2002/95/EC as amended, UE 10/2011 as amended. EU CONEG & it's following amendments. It is safe use in contact with food stuff, Pharmaceuticals & drinking water

Technical Data Sheet

1.4.5 Product is manufactured without the use of Raw Materials of animal origin. Therefore the BSE/TSE issue does not concern this product to the best of our Knowledge. No materials containing latex or natural rubber are used in the Manufacturing, handling & packaging process for this product.

1.5 Color: Transparent (Natural color)

1.6 Quality control:

Product are manufactured, inspected and found to be in compliance with product & quality specification requirements as documented in ISO 9001:2015 & 13485:2016 Quality Management System. Testing of various characteristics in accordance with the valid specification, inspection reviewed & signed for product release.

1.7 Other information

1.7.1 Petri Dishes are suitable for all bacteriological applications.

1.7.2 Vented dishes have an elevated lid that allows for superior gas exchange and evaporation.

1.7.3 Petri dishes retain their shape up to 55°C.

1.7.4 Easy to Stack ,12° inclination 10 stacking filed with 20 ml media from 90mm.

1.7.5 Recommended for single use only

1.8 Basic Features

1.8.1 Temperature range: Stable up to 60°C.

1.8.2 Autoclavability: NO

1.8.3 Chemical resistance: See website (www.abdoslifesciences.com)

1.8.4 Manufactured in Clean room ISO Class 8 (100,000)

1.8.5 Intended use: Widely used in scientific studies of molecular biology, clinical chemistry, biochemistry and any type of scientific laboratory.

Instruction for Use

Product can be used for food contact also.

Storage temp. Max 30°C

Before using the product check crackness, pouch punctured etc.

Check all the information printed on pouches or boxes.

Always wear gloves while handling specimens.

WARNING

DO NOT store strong oxidizing agents.

DO NOT place any plastic labware in a flame.

DO NOT mix any chemicals that may result in a thermic reaction, which can cause product failure

Consult local fire codes prior to storage of flammable liquids.

Technical Data Sheet

UKRP
Obelis UK LTD
Sundford Gate,
East Point Business Park,
Oxford OX4 6LB
United Kingdom
Tel: +44 1491 378012
info@obelis.co.uk

Obelis s.a.
Bd General Wals 53
1030 Brussels Belgium
Tel: +(32) 2 732-59-54
Fax: +(32) 2 732-60-03
mail@obelis.net

NON STERILE
30°C
CE
IVD
EC REG
Obelis s.a.
Bd General Wals 53
1030 Brussels Belgium
Tel: +(32) 2 732-59-54
Fax: +(32) 2 732-60-03
mail@obelis.net

Manufactured & Marketed By :
Abdos Labtech Private Limited
Plot No. 2, Lakshewari Bhagwanpur Roorkee 247667 (UK) India
For Complaints/Feedback Contact Customer Care at above address or
Ph.: +91-011-41081897 or Email us at labtech@abdosindia.com
Lic. No.: MFG/IVD/2024/000158 MADE IN INDIA

Electronic IFU
Scan to visit website

Quality Manager

Abdos Labtech Private Limited

Employee Code 30110

