



Specification

1.07381.0000 Lead(II) carbonate for analysis EMSURE® ACS

Spec. Values		
Assay (complexometric)	≥ 99.0	%
Substances insoluble in dilute acetic acid	≤ 0.02	%
Chloride (Cl)	≤ 0.002	%
Nitrite and Nitrate (as NO ₃)	≤ 0.005	%
Total nitrogen (N)	≤ 0.02	%
Ca (Calcium)	≤ 0.01	%
Cd (Cadmium)	≤ 0.002	%
Cu (Copper)	≤ 0.002	%
Fe (Iron)	≤ 0.002	%
K (Potassium)	≤ 0.02	%
Na (Sodium)	≤ 0.05	%
Zn (Zinc)	≤ 0.003	%

corresponds to ACS

Kristin Haufschildt

Responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature.