

FICHE TECHNIQUE

GRADE ETC

DESIGNATION

- Cartouche d'extraction en cellulose
- Matière : 100% cellulose
- Liant : aucun

APPLICATIONS

- Détermination du taux de lipides dans les produits alimentaires
- Détermination du taux de liant dans le bitume et l'asphalte
- Analyse de pesticides dans les sols
- Extraction de plastifiants dans le PVC
- Extraction de combinaisons organiques à partir des mélanges réactionnels, aliments, substances naturelles, laques et analyse de liant dans les couleurs de peinture, détermination de la poussière et du goudron dans les gaz

SPECIFICATIONS TECHNIQUES

CARACTERISTIQUES	VALEURS	UNITES
Epaisseur	1	mm
Porosité	5-10	µm

DONNEES COMMERCIALES

- Dimensions commerciales de base communiquées pour le diamètre intérieur uniquement

TECHNICAL DATA

GRADE ETC

DESCRIPTION

- Cellulose extraction thimbles
- Material: 100% cellulose
- Binder: no

APPLICATIONS

- Determination of the rate of lipid in food products
- Determination of the rate of binder in bitumen and asphalt
- Analysis of pesticides in soils
- Extraction of plasticizers in PVC
- Extraction of organic combinations out of the reaction mixtures, food, natural substances, lacquers and binder analysis in paint colours, dust and tar determination in gases

TECHNICAL SPECIFICATIONS

CHARACTERISTICS	VALUES	UNITS
Thickness	1	mm
Pore size	5-10	µm

SALES INFORMATION

- Standard sales information given for inner diameter only