

TECHNICAL DATASHEET

Article no.: 7695933

LABSOLUTE[®] pipette tips 500 – 5000 µl


500 pieces in bulk pack

Parameter	Specification
Material	Polypropylene
Colour	transparent
Sterile	NO
Maximum volume	5000 µl
Length	124 mm
Graduation marks	1000, 1500, 2000, 2500, 3000, 3500, 4000, 4500, 5000 µl
Fit	Gilson and Rainin
Certified free from	RNase, DNase and pyrogens

This document has been generated automatically and is valid without signature.

TECHNICAL DRAWING

Article no.: 7695933



This document has been generated automatically and is valid without signature.