


LABSOLUTE® ERLLENMEYER FLASKS

wide neck, with beaded rim

Properties:

- Made of borosilicate glass 3.3
- According to DIN EN ISO 24450
- With beaded rim and labelling field
- Durable, white volume scale
- Dishwasher proof and autoclaveable at 121 °C
- Very high resistance against a wide range of chemicals
- High temperature resistance up to 500 °C

Technical drawing / picture:



The technical drawing is just an example. Technical drawings of all sizes are available on request.

Th. Geyer GmbH & Co. KG

Value table:

Item no.	V ml	Scale	D ml	Ø mm	h mm
7.690 036	25	15 – 25	10	42	75
7.690 037	50	25 – 50	25	51	85
7.690 038	100	50 – 75 – 100	25	64	105
7.690 039	200	100 – 150 – 200	50	79	131
7.690 040	250	100 – 150 – 200 – 250	50	85	140
7.690 041	300	100 – 150 – 200 – 250 – 300	50	87	156
7.690 042	500	200 – 300 – 400 – 500	100	105	175
7.690 043	1,000	400 – 600 – 800 – 1000	200	131	220
7.690 044	2,000	1000 – 1500 – 2000	500	153	280

Description of the abbreviations in the value table:

Item no.	Item number
V	Nominal volume of the Erlenmeyer flask in milliliter (ml)
Scale	Scaling that is printed on the Erlenmeyer flask
D	Division of the scale in milliliter (ml)
Ø	Max. diameter of the Erlenmeyer flask at its widest part in millimeter (mm)
h	Max. height of the Erlenmeyer flask in millimeter (mm)

Further information is part of the declaration of conformity.

Revision 1.0, Stand: 20.07.2017

Th. Geyer GmbH & Co. KG