



Product specification

Product number 2.04737
Packaging number 2.04737.991
Product name Guaiacol Detection Kit

Version 2.1
Date of version 14.06.2018
Page 1/1

Information about composition and property

Appearance and consistency

solid

Odour and taste

unfit for consumption, therefore not applicable

Descriptions

Application

Enzyme assay for detection of beverage spoiling and off-flavour producing Alicyclobacillus spp. in the beverage and fruit compound industry.

Description

40 test tubes each with 2 ml BAT-Induction broth
1 vial with 50 ml Buffer solution
1 vial with 900 µl Hydrogen peroxide Solution (1,3%)
2 vials each with 900 µl Guaiacol-Standard-Solution (1050ppm)
1 vial with Horseradish Peroxidase lyophilised (~90 U)

Other information

Sterilized and ready for use.

Storage and packaging information

Packaging and content: DMD Gaujakol Detection kit, (40 pieces)

Shelf life in the unopened original packaging

At 4-8 °C : 360 day(s)

Storage and handling

Store dry.

Keep away from sources of ignition.

Microbiology

No data available.

Dangerous good and hazardous material information

Disclaimer

This product specification describes the composition and respectively the properties of our product to the best of our knowledge. It replaces all previous versions and is valid without signature.