

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Potassium hexacyanoferrate(III) – ReagentPlus®, ~99%

Product Number:

P8131

CAS Number:

13746-66-2

MDL:


MFCD00011392

Formula:

C₆FeK₃N₆

Formula Weight:

329.24 g/mol



TEST

Specification

Appearance (Color)

Orange to Very Dark Orange and Red to
Very Dark Red and Orange-Red and Red-Orange

Appearance (Form)

Powder or Crystals

Titration by Na₂S₂O₃

≥ 98.5 %

Recommended Retest Period

3 Years

Specification: PRD.3.ZQ5.10000003902

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.