

ICF 400 SUPER- forced air Incubator, capacity 400LT, with touch screen

Item Code: 41100932



DATA SHEET Features

Power supply	230 V
Capacity	400 L
Maximum shelf load	20 kg
Security class	3.1
Dimension	1285 x 800 x 1015 mm
Internal dimensions	1000 x 800 x 500 mm
Display	Touch screen
Minimum distance between shelves	50 mm
Frequency	50 - 60 Hz
Warranty	2 anni
Image 2	https://app.onpage.it/api/storage/6ed595d39e2b3ae3506329a54d8efa69d86f
Working method	Forced ventilation
Name	ICF 400
Regulations	2006/95/EC -2004/108/EC - EN61326-1:2006 - EN61010-1:2010 - EN61010-2-010:2003
Number of shelves	3
Temperature homogeneity at 37 °C	±0,5 °C
Weight	160 kg
Absorbed power	1500 W
Programs	8 programmes of 8 steps each and single step programme
Temperature resolution	0.1 °C
Data transfer	YES, via USB stick
Overtemperature security	Electronic and electromechanical
Working temperature	+5 Ta...80 °C
Heating time (37°C)	55 m
Timer	99 : 59 hh:mm
Ventilation	At 3 speeds (50 %, 75 %, 100 %)
Temperature variation at 37 °C	±0,3 °C
Useful volume	400 L
