

CL00.0619

INDEX

Phenolphthalein indicator

Fenolftaleine indicator
Fenolftaleína indicador
Phénolphtaléine indicateur
Phenolphthalein Indikator
Fenolftaleína indicatore

Titrimetric indicator (pH 8.2-9.8 - colourless to red-violet)
C20H14O4

Mol weight: 318.33 g/mol
Density: 1,30 g/cm³
HS Nr: 29322010
CAS Nr: 77-09-8
UN:
ADR:



Specifications:

Absorption maximum	: 551 - 554 nm
Transition range	: pH 8.2 - 9.8 (colorless - red-violet)
Transition range	: Passes test ACS
Loss on drying	< 1.0%