

UV Lamp

For the detection of fluorescent substances.

The UV lamp can be used for the detection of *Escherichia coli* by means of MUG cleavage in a culture medium containing MUG (Fluorocult® or Chromocult® media, page ...) or in conjunction with Bactident® *E. coli* (Merck, Cat. No. 1.13303.).

Technical Data

Wattage:	4 Watts
Wavelength:	366 nm
Batteries:	5 x 1.5 V baby cells
Weight:	Approximately 400 g
Dimensions:	16 x 9 x 2.5 cm

Ordering Information

Product	Merck Cat. No.	Pack size
UV Lamp	1.13203.0001	