



Automatic Sanitization Module (A.S.M.)

A sound investment for stored water purity

General Description

The ASM (Automatic Sanitization Module) is a device that illuminates the inside of our storage tanks using a bactericidal UV lamp, selected for its germicidal effectiveness.

Key Benefits

- Automatic UV exposure reduces bacteria levels
- Eliminates or limits significantly the need for chemical sanitization
- Optional water detection feature protects lab against water leakages

Ordering Information

| Description | Cat. No. |
|---|-----------|
| Milli-Q® Integral; Milli-Q® Direct, Elix® Advantage, and Elix® Reference water purification systems | TANKASMIN |
| ASM for RiOs™/Elix®/AFS® Essential water purification systems | TANKASMES |



To place an order or receive technical assistance

In Europe, please call Customer Service:

France: 0825 045 645
Germany: 069 86798021
Italy: 848 845 645
Spain: 901 516 645 Option 1
Switzerland: 0848 645 645
United Kingdom: 0870 900 4645

For other countries across Europe, please call: +44 (0) 115 943 0840

Or visit: merckmillipore.com/offices

For Technical Service visit:
merckmillipore.com/techservice

merckgroup.com/life-science

The life science business of Merck operates as MilliporeSigma in the U.S. and Canada.

© 2018 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved. Merck, Milli-Q, Elix, AFS, RiOs, and the vibrant M are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

MK_DS1253EN
2017 - 08107
1/2018