

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com Email USA: techserv@sial.com Outside USA: eurtechserv@sial.com

Product Specification

Product Name:

Copper(II) oxide - nanopowder, <50 nm particle size (TEM)

 Product Number:
 544868

 CAS Number:
 1317-38-0

 MDL:
 MFCD00010979

Formula: CuO

Formula Weight: 79.55 g/mol

CuO

TEST Specification

Appearance (Color)

Appearance (Form)

Complexometric EDTA

Black

Powder

77.0 - 82.6 %

% Copper

X-Ray Diffraction Conforms to Structure

Size < 50 nm

Average Particle Size (by TEM)

Specification: PRD.0.ZQ5.10000096574