

FICHE TECHNIQUE

REF PH114

DESIGNATION

- Papier indicateur
- pH 1.0-14.0

CARACTERISTIQUES

- Support papier imprégné à base de réactifs
- Aspect : Solide (papier)
- Couleur : Jaune /orangé

CONSERVATION

- Au sec et à l'abri de la lumière

APPLICATIONS

- Analyses de solutions en milieu faiblement acide, neutre ou légèrement alcalin

PRECISION DE LA MESURE

- La mesure est étalonnée par saut de 1.0 unité de pH
- Le papier réactif réagit par méthode colorimétrique ; le résultat obtenu doit être comparé à l'échelle fournie avec chaque article.

DONNEES COMMERCIALES

- Disponible en rouleaux, en bandelettes et en livrets

TECHNICAL DATA

REF PH114

DESCRIPTION

- Indicator paper
- pH 1.0-14.0

CHARACTERISTICS

- Paper support impregnated with reagents
- Appearance: Solid (paper)
- Colour: Yellow / orange

STORAGE

- Dry and light-free storage

USE

- Analysis of solutions in weakly acidic, neutral or slightly alkaline medium

ACCURACY OF MEASUREMENT

- The measurement is calibrated every 1.0 pH unit.
- The reagent paper reacts by colorimetric method; the result obtained must be compared to the colorimetric scale given with each article.

SALES INFORMATION

- Available in rolls, strips and booklets