

FICHE TECHNIQUE

REF PH590

DESIGNATION

- Papier indicateur
- pH 5.0-9.0

CARACTERISTIQUES

- Support papier imprégné à base de réactifs
- Aspect : Solide (papier)
- Couleur : Jaune

CONSERVATION

- Au sec et à l'abri de la lumière

UTILISATIONS

- Analyses de solutions en milieu faiblement acide, neutre ou légèrement alcalin

PRECISION DE LA MESURE

- La mesure est étalonnée par saut de 0.5 unité de pH
- Le papier réactif réagit par méthode colorimétrique ; le résultat obtenu doit être comparé à l'échelle fournie avec chaque article.

DONNEES COMMERCIALES

- Disponible en rouleaux, en bandelettes et en livrets

TECHNICAL DATA

REF PH590

DESCRIPTION

- Indicator paper
- pH 5.0-9.0

CHARACTERISTICS

- Reagent-impregnated paper backing
- Appearance: Solid (paper)
- Colour: Yellow

STORAGE

- Dry and light-free storage

APPLICATIONS

- Analysis of solutions in weakly acidic, neutral or slightly alkaline medium

PRECISION DE LA MESURE

- The measurement is calibrated every 0.5 pH unit.
- The reagent paper reacts by colorimetric method; the result obtained must be compared to the colorimetric scale given with each article.

SALES INFORMATION

- Available in rolls, strips and sheets