

BagPipet® & BagTips®

Fixed-volume pipet & straw tips

Ref 251 091 (BagPipet®) - Ref 252 019 (BagTips® Regular) - Ref 252 024 (BagTips® Jumbo)



Standards:



Manufactured under:



Product made for INTERSCIENCE by Interlab, an ISO 9001 certified company



Distribué par :

Z.A de Gesvrine - 4 rue Képler - B.P.4125
44241 La Chapelle-sur-Erdre Cedex - France
t. : +33 (0)2 40 93 53 53 | f. : +33 (0)2 40 93 41 00
commercial@humeau.com

BagPipet® is an autoclavable pipet with fixed volume (0.1 mL / 0.9 mL / 1 mL). Used with BagTips® a single use rigid straw tip, available in 2 sizes (Regular 19 cm / Jumbo 24 cm), BagPipet® intakes any product, even very viscous products: milk, cream, yoghurt, cosmetics...The long BagTips® reaches the bottom of any bag.

TECHNICAL SPECIFICATIONS

- Product made for INTERSCIENCE by Interlab, an ISO 9001 certified company
- Designed and made in France

BagPipet®

Ref: 251 091

- Fixed volume pipet
- 3 settings: 0.1 mL / 0.9 mL / 1 mL
- Calibrated according to the GLP "Good Laboratory Practice"
- Autoclavable
- Ejection of disposable straw
- Easy calibration
- Light and ergonomic
- Also available: variable volume button from 10 to 1000 µL (Ref 251 080)
- Box (w x d x h): 30 x 10 x 6 cm, weight: 250 g
- 1 year warranty
- In compliance with: ISO 7218, ISO 6887-1, ISO 4833-2 and FDA BAM (Bacteriological Analytical Manual)

DELIVERED WITH

- 1 user's manual
- 1 pouch BagTips® Regular
- 1 fixed volume button 1000 µL (on the pipet)
- 1 holder for key
- 2 fixed volume buttons (900 µL, 100 µL)
- 2 O-Rings
- 1 Allen wrench
- Oil

BagTips® Regular

Ref: 252 019

- Nominal volume: 1.1 mL
- Treated by Gamma ray 5 to 12 kGy, with certificate
- Length: 190 mm
- Quantity per box: 1000
- Box (w x d x h): 27.5 x 23 x 11 cm, weight: 1.2 kg

BagTips® Jumbo

Ref: 252 024

- Nominal volume: 2.2 mL
- Treated by Gamma ray 5 to 12 kGy, with certificate
- Length: 240 mm
- Quantity per box: 1000
- Box (w x d x h): 27.5 x 23 x 11 cm, weight: 1.5 kg