

Material 20839.297  
 Material description 2-Propanol  
 Grade GPR RECTAPUR

CAS Number 67-63-0  
 Molecular formula  $(\text{CH}_3)_2\text{CHOH}$   
 Molecular mass 60.10

| Characteristics                | Specifications |
|--------------------------------|----------------|
| Assay                          | ≥ 99.0 %       |
| Acidity                        | Passes test    |
| Acidity or alkalinity          | Passes test    |
| Benzene and related substances | Passes test    |
| IR Spectrum                    | Passes test    |
| Peroxides                      | Passes test    |
| Boiling point                  | 81.5 - 82.6 °C |
| Density (20/20)                | 0.785 - 0.789  |
| Density (25/25)                | 0.783 - 0.790  |
| Refractive index (20°C)        | 1.376 - 1.378  |
| Evaporation residue            | ≤ 20 ppm       |
| Water                          | ≤ 0.5 %        |

## Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

Isabelle Habay, Head of Laboratory - Briare  
 VWR International S.A.S.; Z.I. de Vaugereau; FR-45250  
 Briare; France

For Professional use in Laboratory or Manufacturing. Not for use as an Active Pharmaceutical Ingredient or Food or Animal Feed. Suitability and intended use of the product remains the responsibility of the user.