


CATALOG NO.	5245-0040	5245-0050			
SIZE CODE	D	E			
MAX. PIPET LENGTH (MM/IN.):	406/16	610/24			
MATERIAL	BASE	BLUE HOPE			
	BODY	HOPE			
	TUBING *	PVC			
		INCHES	MM	INCHES	MM
A (BASE DIA.)	9.0	229	9.0	229	
B (OVERALL WIDTH)	10.7	272	10.5	267	
D (NOM. OVERALL HEIGHT)	30	762	38	965	
H (NOM. HEIGHT)	27	686	35	889	
NOM. WT. (gms)	1800		1700		

* INCLUDES (2) PIECES OF TUBING AS FOLLOWS:
(1) 1" I.D. x 1 1/4" O.D. x 18" LONG
(1) 1/2" I.D. x 3/4" O.D. x 48" LONG

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Ponorono Creek Drive
Rochester, NY 14625-2385, USA

MADE	MCC	APP'D	PVC	PSS 00-5245-00X0	REV	0995
EX'D	PVC		9/14/95			

DESCRIPTION:

PIPET WASHER-RINSER