

27. Jul 2015



## Sterilization indicator tape

### Steril.indicator tape 50 m x19 mm crepe

Minimum order quantity / quantity interval: 1 roll

### Product description

Self-adhesive crêpe paper containing thermo-active pigments. The white colored strips turn brown after autoclaving (121 °C: after 20 min. approx., 134 °C: after 5 min. approx.).

### Specification

length	m	50
width	mm	19

Distribué par :

Z.A de Gesvine – 4 rue Képler – B.P.4125  
 44241 La Chapelle-sur-Erdre Cedex – France  
 t. : +33 (0)2 40 93 53 53 | f. : +33 (0)2 40 93 41 00  
 commercial@humeau.com



w w w . h u m e a u . c o m