



Produit	31627.290
Description du produit	Sodium hydroxyde 1 mol/l (1 N) en solution aqueuse
Grade / qualité	AVS TITRINORM Reagent Ph.Eur. chapter 4.2.2, Reagent USP-NF
Numéro CAS	1310-73-2
Formule moléculaire	NaOH
Masse moléculaire	40.00

Information additionnelle

Analyses	Spécifications
Titre(20°C;valeur vraie 0,2% exactitude)	0.998 - 1.002 Mol/l
Aspect	Liquide limpide incolore
Conforme à Réact. Ph.Eur.	Conforme
Conforme à Réact. USP-NF	Conforme
Traçable selon NIST	Confirmé

Signature	Information additionnelle
Nous certifions que ce lot est conforme aux spécifications listées ci-dessus. Ce document a été produit électroniquement et est validé sans signature. Isabelle Habay, Head of Laboratory - Briare VWR International S.A.S.; Z.I. de Vaugereau; FR-45250 Briare; France	La valeur garantie est rattachée à un matériau de référence primaire (SRM) de chez National Institute of Standards and Technology (NIST).

For Professional use in Laboratory or Manufacturing. Not for use as an Active Pharmaceutical Ingredient or Food or Animal Feed. Suitability and intended use of the product remains the responsibility of the user.