

---

# XS Basic pH 4.01 / 25°C, 55ml bottle (red) Test solution



## TECHNICAL DATA SHEET

### Features

---

|               |                                              |
|---------------|----------------------------------------------|
| Accuracy      | ± 0.02                                       |
| Colour        | Pink                                         |
| Composition   | Potassium phthalate acid, potassium chloride |
| Package       | Disposable bottle 55 ml                      |
| Description   | ±0.02 / 25°C - Red                           |
| Name          | pH 4.01                                      |
| Nominal value | pH 4.01 @ 25 °C                              |
| Volume        | 55 ml                                        |
| hazards       | Non-hazardous. No ADR restrictions           |

---