
XS Basic pH 7.00 / 25°C, 55ml bottle (green) Test solution

Code: 32383133



TECHNICAL DATA SHEET

Features

Accuracy	± 0.02
Colour	Green
Composition	Monobasic potassium phosphate (CAS 7778-77-0), sodium dibasic phosphate dihydrate (CAS 10028-24-7), green colouring
Package	Disposable bottle 55 ml
Description	±0.02 / 25°C- Green
Name	pH 7.00
Nominal value	pH 7.00 @ 25 °C
Volume	55 ml
hazards	Non-hazardous. No ADR restrictions
