
1X250ml XS Solution pH 7.00 ±0.01/25°C (green), with certificate traceable to NIST

Item Code: 51100043



DATA SHEET

Features

Accuracy	± 0,02
Certificate	Yes
Colour	Green
Composition	Monobasic potassium phosphate (CAS 7778-77-0), sodium dibasic phosphate dihydrate (CAS 10028-24-7), green colouring
Pack	250 ml with integrated calibration compartment
Description	±0.02 / 25°C - Green
Name	pH 7.00
Nominal value	pH 7.00 @ 25 °C
Volume	250 ml
