



Thermo Scientific™ *Brilliance*™ MRSA

User Guide

presumptive identification of Methicillin-resistant *Staphylococcus aureus* (MRSA)



Inoculate plate. Incubate at 37°C for 18-24 hrs.

Read results.



Positive



Blue

Negative



Pink, purple or white colonies

Identifications are presumptive and should be confirmed. Please note, organisms with an atypical enzyme pattern may give anomalous reactions on Thermo Scientific™ *Brilliance*™ MRSA Agar.

thermoscientific.com/preparedmedia

© 2014 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries.

Thermo
SCIENTIFIC

A Thermo Fisher Scientific Brand