



# Specification

1.09204.0100 Giemsa's azur eosin methylene blue solution for microscopy

Specification	
Suitability for microscopy	passes test
Hematological staining	
Specification	
Erythrocytes	reddish-brownish
Nuclei	violet
Eosinophilic granules	reddish to red-brown
Neutrophilic granules	light violet
Lymphocyte cytoplasm	blue
Monocyte plasma	grey-blue
Thrombocytes	violet
Basophilic granules	dark violet
Blood parasite nuclei	light red
Histological staining	
Specification	
Nuclei (Tissue section)	blue to dark blue
Cells	blue to dark blue
Collagen	pale blue
Osteoid	pale blue
Eosinophilic granula (Tissue section)	red
Acidic mucopolysaccharides	reddish-violet
Mast cell granules	reddish-violet
Cartilage matrix	reddish-violet
Acidophilic materials	orange-red
Helicobacter pylori	blue to dark blue

Dr. Timm Schlegelmilch  
Responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature.