



**SUPPLEMENTI SELETTIVI ED ARRICCHIMENTI  
SELECTIVE SUPPLEMENTS AND ENRICHMENTS**

## **ROSOLIC ACID SUPPLEMENT**

Supplemento per l'isolamento di coliformi fecali.  
*Supplement for faecal coliforms isolation.*

### **MODALITA' D'USO**

Ricostruire asetticamente il contenuto di una fiala con 5 ml di acqua distillata . Agitare fino a completa dissoluzione ed aggiungere a 500 ml di m-FAECAL COLIFORM Agar (m-FC) (Codice 611012-6110125 - 621012). Miscelare bene. Portare ad ebollizione. Non sterilizzare. Raffreddare a 50°C. Distribuire in piastre petri.

### **INSTRUCTION FOR USE**

*Aseptically reconstitute the content of one vial with 5 ml of sterile distilled water . Mix to complete dissolution and add to 500 ml of m-FAECAL COLIFORM Agar (m-FC) (Code 611012 – 6110125 - 621012). Mix well. Bring to boil. Do not sterilize. Cool to 50°C. Distribute into petri dishes.*

### **FORMULA (grammi per fiala)**

Acido Rosolico..... 50 mg

### **FORMULA (grams per vial)**

Rosolic Acid..... 50 mg

### **CONTENUTO**

10 Fiale liofili di ROSOLIC ACID Supplement

### **CONTENT**

10 Lyophilic vials of ROSOLIC ACID Supplement

### **CONSERVAZIONE**

2-8°C

### **STORE AT**

5-12°C

### **VALIDITA'**

3 anni

### **SHELF LIFE**

3 years

<b>PRODOTTO / PRODUCT</b>	<b>CODICE / CODE</b>	<b>CONFEZIONE / PACKAGING</b>
ROSOLIC ACID SUPPLEMENT	81029	10 Fiale per 500 ml di terreno/ 10 Vials for 500 ml of medium
m-FAECAL COLIFORM AGAR (m-FC)	611012	500 g
m-FAECAL COLIFORM AGAR (m-FC)	6110125	5 K g
m-FAECAL COLIFORM AGAR (m-FC)	621012	100 g

### **BIBLIOGRAPHY**

1. APHA (1985). *Standard Methods for the Examination of Water and Wastewater*, 16<sup>th</sup> Edition.
2. Gildreich, Clarck, Huff and Bert (1965). *J. Am. Waters Works Ass.* 57, 208.

Rev.1/01.02.2007