

PRODUCT SPECIFICATION

Specification	Product Code	BNP05A
	Description	5L Aspirator, Heavy Duty
	Material	Polypropylene
	Batch No.	N/A
	Colour	Opaque
	Autoclavable	Yes
	Pack quantity	1

Material Information	Temperature Range	-20 to+135 ⁰ C
	Autoclave Temperature	121 ⁰ C
	Chemical Resistance	Good to Excellent
	Further details:	https://www.dwkLtd.com/en/chemicalcompatibility

The above product is compliant to the following standards only if marked 'Yes'

Compliance	Food safe (EU/FDA)	YES- EU 10/2011 / FDA CFR 21 177
	BSE / TSE free	YES- No tallow from animal sources are used in the manufacture of the finished product.
	REACH Compliant	YES- Compliant according to 1907/2006
	LATEX	YES- No Materials containing latex natural rubber are used in the manufacturing, handling, and packaging process.

DWK Life Sciences Limited hereby confirms that the product specified above is inspected according to all relevant specifications and procedures laid out in our Quality System, in accordance with our registration to BS EN ISO9001 and BS EN ISO13485.

Signed on behalf of DWK Life Sciences Limited by:

Signed: M.Ashley
Quality Technician
Date: 12 Oct 2020