



Product Specifications

CL00.0159

Ascorbic acid (L(+)) v.p.

Ascorbinezuur (L(+)) z.z.
Acido ascorbico (L(+)) e.p.
Acide ascorbique (L(+)) t.p.
Ascorbinsäure (L(+)) Rst.
Acido ascórbico (L(+)) t.p.

For laboratory use
99+% C₆H₈O₆

Mol weight: 176.13 g/mol
Density: 1,65 g/cm³
HS Nr: 29362700
CAS Nr: 50-81-7
UN:
ADR:

Specifications:

Assay	>99%
Chloride (Cl)	<0.005%
Sulfate (SO ₄)	<0.002%
Oxalic acid	<0.2%
Methanol	<0.1%
Arsenic (As)	<0.0003%
Copper (Cu)	<0.0005%
Iron (Fe)	<0.0002%
Mercury (Hg)	<0.0001%
Lead (Pb)	<0.0005%