



Product Specifications

CL05.2615

TITREX

Sulfuric acid 0.5 mol/l

Acido sulfúrico 0.5 mol/l
Zwavelzuur 0.5 mol/l
Acide sulfurique 0.5 mol/l
Schwefelsäure 0.5 mol/l
Acido solforico 0.5 mol/l

For laboratory use, Titrations

49.04 g H₂SO₄ / l H₂O = 1 N (±0.002/20°C)

Traceability: Reference standards acc. NIST SRM 84I

Mol weight: 98.08 g/mol

Density: 1,03 g/ml

HS Nr: 28070000

CAS Nr: 7664-93-9

UN: 3264

ADR: 8,III



Specifications:

CHEM LAB NV

Industriezone "De Arend " 2 B-8210 ZEDELGEM - BELGIUM

Tel.: +32 50 28 83 20 Fax. +32 50 78 26 54 e-mail : info@chem-lab.be Web : www.chem-lab.be