

Specification

Aquabator-Clean™ (100X)

A9390

Product Code:	A9390
Product Name:	Aquabator-Clean™ (100X)
Headline Comment:	• Aquabator-Clean™ is a registered trademark of AppliChem GmbH.
Specifications:	: • Disinfectant solution for ordinary water baths (not for CO₂-incubators!) : • for the prevention of microbial growth in water baths : • 100X concentrated solution (Use 10 ml per 1 liter of water.)
Hazard pictograms	
WGK:	2
Storage:	2 - 8°C protected from light
Signal Word:	Attention
GHS Symbols:	GHS07
H Phrases:	H319
P Phrases:	P264 P280 P305+P351+P338 P337+P313
CS:	38220000

AppliChem GmbH

Ottoweg 4 • D-64291 Darmstadt • Phone +49 6151 9357 0 • Fax +49 6151 9357 11 • info.de@itwreagents.com • www.itwreagents.com
 CEO Joan Roget • Commerzbank Darmstadt • Bank 508 800 50 • Account 0186989900 IBAN DE24 5088 0050 0186 9899 00 • Swiftcode DRESDEFF508 • Finanzamt Darmstadt 07 228 16476 • Register court Darmstadt HRB Nr. 7340