


BL 303



- **Backlight LCD display**
- **Rechargeable battery**
- **RS 232 port**
- **Capacity and stability indicator**
- **21 different measures units**
- **Supplied with external check weight**
- **Supplied with calibration report .**
(only BL 303)



Keypad and display

BL 602



BL 2002



BL 3001 / BL 6001



Code	Model	Capacity	Division	Reproducibility	Linearity	Calibration	Plate dimension
042.000000.03	BL 303	300 g	0,001 g	0,001 g	± 0,003 g	External	Ø 118 mm
042.000000.06	BL 602	600 g	0,01 g	0,01 g	± 0,03 g	External	Ø 118 mm
042.000000.04	BL 2002	2000 g	0,01 g	0,01 g	± 0,02 g	External	188 x 148 mm
042.000000.08	BL 3001	3000 g	0,1 g	0,1 g	± 0,2 g	External	188 x 148 mm
042.000000.07	BL 6001	6000 g	0,1 g	0,1 g	± 0,2 g	External	188 x 148 mm