



designed for scientists



C 5010 Decomposition vessel, standard

/// Data Sheet

Material: stainless steel

With fix ignition wire

Can be converted to use the C 14 Combustible crucible



designed for scientists

Technical Data

Useable volume max. [ml]	260
Pressure max. [bar]	230
Test pressure [bar]	330
Diameter [mm]	57
Working temperature max. [°C]	50
Oxygen pressure max. [bar]	40
Amount of energy max. [J]	40000
Material in contact with medium	stainless steel 1.4462
Height [mm]	136
Weight [kg]	2.1

