

# BM 113

## Moisture Analyzer



### Data sheet

Capacity	110 g
Division	0,001 g
Repeatability	0,15 %
Min. weigh	0,008 g
Min. weight for test	1 g
Time range	1-99 min – 1 min step
Temperature range	40 – 199 °C
Heating control	Timing – Auto - Manual
Heating control	Standard – Fast - Soft
Calibration	External
Plate dimension	Ø 90 mm
Lamp	Halogen Lamp 400W
Memory	15 programs
Communication	Serial Port RS-232
Dimensions (mm)	340 * 200 * 180
Net weight	9 Kg
Power supply	220 V – 50Hz

Distribué par :

Z.A de Gesvrine – 4 rue Képler – B.P.4125  
44241 La Chapelle-sur-Erdre Cedex – France  
t. : +33 (0)2 40 93 53 53 | f. : +33 (0)2 40 93 41 00  
commercial@humeau.com

